



Lifestyle P2

Portable
Oxygen
Concentrator



Specifications

Model:	P2
Size:	8.7" x 3.4" x 6.3"
Weight:	4.37 lbs
Warm up time:	2 minutes
Oxygen Flow:	Pulse Oxygen Supply
Oxygen Concentration:	87%-96%
Oxygen Flow:	Pulse dose delivery system, 5 flow settings
AC Power Supply:	100-240v, 50-60hz (auto sensing to allow for worldwide use)
Rechargeable Battery:	14.4VDC
Battery Rated Capacity:	6.8Ah
Display Screen:	2.8"
Battery:	6800mAh, duration up to 4 hours Approx. 4 hours charge time with AC
Operation:	Simple control functions and easy-to-read LCD display
Use:	Designed for 24/7 use
Warranty:	Three years, excluding battery and sieve bed


Lifestyle Mobility Aids | 3200 Tyrone Blvd N | St. Petersburg, FL 33710
Email: contactus@lifestylemobilityaids.com | Website: www.lifestylemobilityaids.com


Compact & Lightweight





Powerful Compressor




- 

High Efficiency
- 


Low Noise
- 


Excellent Heat Dissipation
- 


Long Life





Wherever you go Oxygen follows


- 

4.4 lbs
Lightweight
- 

Equivalent
to 5 Liters
- 

High Absorption
Sieve Bed
- 

Rechargeable
Battery
Duration up to 4 hrs
- 

93% O₂
Concentration
- 

Small and
Convenient



External Battery Charger*



Conveniently charge a backup battery
while the device is in use.

Charging Time: 2.5 hours
*Additional purchase.



User Control

