



# P2 Portable Oxygen Concentrator



ITEM #	DESCRIPTION
P2	1L Portable Oxygen Concentrator

## FEATURES

- 1 Compact & lightweight.
- 2 Conveniently charge a backup battery while the device is in use. Charging time 2.5 hours. (additional purchase).
- 3 High efficiency, low noise. Excellent heat dissipation.
- 4 93% O<sub>2</sub> concentration.
- 5 High absorption sieve bed.
- 6 Equivalent to 5L.

## SPECIFICATIONS

Size	8.7" x 3.4" x 6.3"
Weight	4.37lbs
Warm up time	2 minutes
Oxygen Flow	0.2 - 1L Pulse Oxygen Supply
Settings	Five Flow Settings
Oxygen Concentrator	87%-96%
AC Power Supply:	100-240v, 50-60hz (auto sensing to allow for worldwide use)
Rechargeable Battery	14.4VDC
Battery Rated Capacity:	6.8Ah
Display Screen	2.8"
Battery	6800mAh, duration up to 4 hours Approx. 4 hours charge time with AC
Operation	Simple control functions and easy-to-read LCD display
Use	Designed for 24/7 use

## WARRANTY

Three years, excluding battery and sieve bed

V1: 12/9/20