



Change Your Oxygen, Change Your Life!

YOUR GUIDE TO RESPIRATORY FREEDOM



Call 844-977-1364 for a free consultation



TAKE A DEEP BREATH. THAT'S WHAT FREEDOM FEELS LIKE.

If you have respiratory issues and require oxygen therapy, get ready for a welcome breath of fresh air. For 20 years, Inogen has helped patients worldwide replace bulky oxygen tanks with lightweight, portable oxygen concentrator systems that deliver superior health benefits with the extra advantage of greater freedom.

From the convenience of our popular Inogen One to innovative new solutions for those who need more specialized care, our mission is to help you receive all the benefits of oxygen therapy without limiting the lifestyle you love. This guide will introduce you to the current Inogen range. To confirm the right option for you, call **844-977-1364** to speak with an oxygen specialist.





WHY WE SELL AND SERVICE INOGEN PORTABLE OXYGEN SYSTEMS

Our story began when Mae, the beloved grandmother of one of our founders, was diagnosed with COPD. Over time, she became reliant on bulky oxygen bottles 24 hours day, making her feel trapped at home. She knew there had to be a better way, and asked her granddaughter why there wasn't a system able to safely deliver oxygen in a portable size. That single question drove us to innovate the entire industry. It inspired us to create a single oxygen concentrator that could be used at home, away, and for travel, the Inogen One.

When Mae tested the Inogen One prototype in 2004, she replied, "This is better than anything I could have ever dreamed of." Since then, Inogen has gone on to become the global leader in portable oxygen therapy, improving the lives of nearly 1 million patients worldwide.

INOGEN ONE.

Portable oxygen to help you enjoy life without limits.


Inogen One is an all-in-one portable oxygen concentrator designed for superior mobility and ease. It replaces lugging around heavy oxygen tanks — or the need for any tanks at all — because it converts the air around you into medical-grade oxygen. All in a system small enough to discreetly carry in a shoulder bag.

While Inogen One is one of the most powerful pulse dose technology systems, it's also super light. Take it anywhere you go with total confidence and freedom. Charging is also a snap, with the ability to plug into a wall outlet, car DC outlet, or on the go with an external battery charger.

Three Inogen One systems are available, with up to six oxygen flow settings. Call us to find the model that's ideal for you.

Freedom Takes Flight

All Inogen One portable systems meet FAA requirements for air travel.



ACCESSORIES TO MAKE YOUR INOGEN ONE EVEN MORE PERSONAL.

As an Inogen One customer, you can select from a wide range of accessories to further enhance comfort, ease, personality and more. Here are just a few of our most popular items.

AC and DC power cables

Additional batteries

External battery chargers

Backpacks



△ △ We always used to travel, and I was

anywhere, even overseas.

afraid those days were over. Thanks to the Inogen One, I'm still able to travel

INOGEN at HOME

WHISPER QUIET OXYGEN THERAPY TO HELP YOU **REST EASY.**

Patients who suffer from chronic respiratory distress often require a continuous flow of oxygen, especially while they sleep at night. Our Inogen At Home system combines powerful, secure 5-liter per minute oxygen delivery with technology that not only removes noise, but greatly reduces size and power consumption.

At just 18 pounds — half the weight of some other oxygen concentrators in use today — Inogen At Home is a compact unit you can easily and safely keep at your bedside. In fact, it's so light you can even use it during the day as a companion around the house. All with incredible energy efficiency to help you save more over time.

Standard Home Oxygen Concentrator

- X Heavy & bulky
- X Very loud
- X High power consumption








- Half the weight!
- ✓ Only 18 pounds
- **Extremely quiet**
- Low power consumption

Call 844-977-1364 to reclaim your freedom



SAVE \$300 OR MORE ON YOUR ELECTRIC BILL EACH YEAR.*

Thanks to its proprietary energy technology, Inogen At Home consumes less electricity than most older, standard home oxygen concentrators. Because it's designed for 24/7 use, the reduced power consumption can significantly add up.




^{*} Calculations are based on published US average hourly kilowatt rates, and assume stationary oxygen concentrator usage of 20 hours per day. A flow setting of 2 liters per minute was used as the Inogen At Home power draw, with the published specifications of a common competitive stationary concentrator used as the standard market concentrator power draw.

freedom. With the Inogen, my quality of life has increased a thousand-fold.

—Claudia H.

Powerful oxygen therapy in a compact size.



^{**} See the Inogen At Home Limited Warranty Statement for full warranty information.

DISCOVER THE RIGHT INOGEN SYSTEM FOR YOU.









Inogen ROVE 6

Inogen One G5

Inogen One G4

Inogen One 🔞



WEIGHT	4.7 lb	4.7 lb	2.8 lb	4.8 lb
OXYGEN FLOW (liters/min)	Pulse-dose, Intelligent Delivery Technology® 1-6	Pulse-dose, Intelligent Delivery Technology® 1-6	Pulse-dose, Intelligent Delivery Technology® 1-3	Pulse-dose, Intelligent Delivery Technology® 1-5
BATTERY LIFE	Single battery: up to 6.5 hours Double battery: up to 13 hours	Single battery: up to 6.5 hours Double battery: up to 13 hours	Single battery: up to 2.7 hours Double battery: up to 5 hours	Single battery: up to 4.7 hours Double battery: up to 10 hours
POWER SENSING	AC Power Supply 100-240 V, 240 W, 50-60 Hz Auto sensing to allow for worldwide use	AC Power Supply 100-240 V, 240 W, 50-60 Hz Auto sensing to allow for worldwide use	AC Power Supply 100-240 V, 240 W, 50-60 Hz Auto sensing to allow for worldwide use	AC Power Supply 100-240V, 50-60Hz Auto-sensing for worldwide use. DC Power for mobile use
WARRANTY	3-Year Warranty + Extended Warranty Available			2-Year Warranty
USE	Day, portable, travel	Day, portable, travel	Day, portable, travel	Day, portable, travel
SIZE	8.5" L x 3.26" W 7.19" High	8.15" L x 3.26" W 7.19" High	5.91" L x 2.68" W 7.2" High	8.75" L x 3" W x 8.25" High
NOISE	38 dBA**	38 dBA**	40 dBA**	39 dBA**



OUR CUSTOMERS ALWAYS HAVE THE FINAL WORD.

- My Inogen system is so much lighter than canisters, and gives me hours of freedom to do what I want. I also love the charging options, and the Inogen Connect app. I only wish I found Inogen sooner.
 - Jerry L.
- As someone with COPD, having Inogen at my side makes me feel safe. Last Christmas, my family took me to Walt Disney World, and I had the best day of my life. Thank you for giving me the ability to live my life again.
 - Situ G.
- 1'm beyond satisfied with my Inogen portable and plug-in concentrators. The portable one lets me keep doing the work I enjoy, and walks and light exercise, while the plug-in is so quiet and lightweight. It travels with me to visit family, and on trips with my husband. I'm so thankful.
 - Dehra F



-Thomas B.

Freedom is the word I think about when someone asks why I chose Inogen. It gives me the freedom to travel ... the freedom to do the things I want to do ... to enjoy life. Inogen is simply the best portable oxygen partner there is.













Discover how people like you have regained their freedom with Inogen. www.oximedical.com



Call 844-977-1364 for a free consultation