

The OxyQ7 - The Latest Innovation in Oxygen Therapy

Meet the OxyQ7

The **OxyQ7** Portable Oxygen Concentrator is designed to give you the freedom to live life on your terms without compromising on performance. Combining advanced oxygen delivery technology with a lightweight, easy-to-carry design, the OxyQ7 delivers reliable support wherever your day takes you. With more pulse flow settings and higher oxygen output than its closest competitors, it provides powerful, responsive therapy you can trust. Whether you're at home, traveling, or enjoying time with loved ones, the OxyQ7 offers dependable performance, greater flexibility, and confidence in every breath.



Mobility Without Limits

The **OxyQ7's** versatile carrying system is built for motion. It comes with a sleek, stylish bag designed to keep up with your lifestyle — whether you're exploring outdoors, navigating the airport, or enjoying a busy day out.

- Wear it comfortably over your shoulder
- Switch to backpack straps for hands-free freedom
- Secure it with the belt strap for a lightweight, low-profile feel

However you move, the OxyQ7 moves with you — delivering the ultimate combination of portability, mobility, and confidence on the go.



Backpack Straps



Belt Strap



Shoulder Strap

The OxyQ7 is Available Exclusively at OxiMedical

Advanced Technology That Sets It Apart

The **OxyQ7** is designed to deliver a higher level of performance compared than other portables in its class!

It provides up to **20% more oxygen** output than its closest competitors.

It also comes with a pulse oximeter — an essential tool for monitoring your oxygen saturation and staying informed about your respiratory health.

In addition, the concentrator itself continuously **monitors your breath rate**, helping alert you to potential changes so you can stay proactive and aware. By combining enhanced oxygen delivery with intelligent monitoring features, the OxyQ7 offers advanced support designed to give you greater confidence in your daily oxygen therapy.



Flexible Power for Every Situation

The **OxyQ7** is designed to keep up with wherever life takes you.

The **home power supply** plugs into any standard AC outlet for easy charging.

The **car power supply** connects to the DC outlet in your car, RV, or even a boat — perfect for travel and long days out.

An **external battery charger** allows you to keep an extra battery charging while you're on the go, so when you return, a fully powered backup is ready to switch out without missing a beat.



Pulse Flow Setting	Battery Life
1	Up to 8 Hours
2	Up to 5.6 Hours
3	Up to 3.4 Hours
4	Up to 2.3 Hours
5	Up to 2 Hours
6	Up to 1.6 Hours
S	Up to 1.4 Hours

Everything You Need. Ready When You Are.

The OxyQ7 is supported by a full range of travel-ready accessories and power options — featuring long lasting batteries, an external battery charger, versatile carrying solutions, and multiple power supply options. However your day unfolds, the OxyQ7 has everything available to keep you powered, prepared, and on the move.



- A. OxyQ7 Portable Oxygen Concentrator
- B. Carrying Case
- C. Accessory Case
- D. Belt Carrying Strap
- E. Shoulder Strap
- F. Backpack Straps
- G. Extra Battery
- H. Car Charger
- I. Wall Charger
- J. External Charger
- K. Pulse Oximeter
- L. Owners Manual
- M. Cannula

OxyQ7 Specs



Pulse Settings
1 - 7



Quiet Operation
39 Decibels



Extensive Battery Life
Up to **8 Hours**



Responsive
Delivery



Home Charger
Car Charger
External Charger



Breath Rate
Detection



Very Lightweight
4 Pounds
With Battery



FAA Approved

Customers are Raving about the OxyQ7



"She loves the weight and LOVES the carry bag — thinks it looked bougie."

"It's super lightweight for the power it puts out. Thanks for keeping me moving."

"I like the versatility of the bag. I tried the backpack and belt options and like wearing it as the backpack the best on my walks."

"Even on a setting of 1 it gave me more oxygen than my tanks."

"She's loving it. It's giving her way more freedom."



The Portable Oxygen Store
OxiMEDICAL

Call Now!

1-888-643-4921

www.OxiMedical.com

We're With You Every Step of the Way.

Purchasing an oxygen concentrator is an important decision — and you shouldn't have to navigate it alone. Our team is here to make the process simple, supportive, and stress-free.

We can check to see if **Medicare can help offset the cost** of your concentrator. Our team handles the paperwork and works to cut through the red tape for you — so you don't have to.

We offer simple approval financing with payments starting as low as **\$59 per month**. Applications can be taken over the phone, and in most cases we can give you instant approval.

Every concentrator purchased from us includes **free lifetime phone support**. Our experienced respiratory specialists are just a call away whenever you need help, answers, or reassurance.

Your OxyQ7 includes **accidental damage protection** for added peace of mind. If your concentrator is dropped or damaged, we'll replace it with a working unit right away — no waiting weeks for repairs or service.

More Than a Supplier — A Partner in Your Care

When you choose OxiMedical, you're not just purchasing equipment. You're gaining a team committed to your mobility, independence, and confidence.